






EDIBON
Technical Teaching Equipment

Computer Controlled Thermal Radiation Unit
TRTC



PROCESS DIAGRAM AND ELEMENTS ALLOCATION





SPECIFICATIONS

Items supplied as standard

① **TRTC. Unit:**

Anodized aluminium structure. Panels and main metallic elements in stainless steel.
Diagram in the front panel with similar distribution that the elements in the real unit.
This unit consist on a metal plate with a resistance at one side and a lamp in another side.
Lengthwise of the metal plate you can place the elements supplied with this unit.
7 Temperature sensors.

This unit is provided with elements for light experiments and radiation experiments:

Light elements:

Lux meter. It has an output to the interface and you can read the value from PC by serial port.
Photodetector.

5 Different grey natural density filters.

Filters porthole.

Radiation elements:

Radiometer:

Planes Surfaces (each one contains a temperature sensor):

Polished Aluminium.


Anodized Aluminium.

Brass.

2 Black Body.

Variable Slit or Aperture, that allows to regulate the area of the radiation.

Wattmeter



TRTC. Unit

② **TRTC/CIB. Control Interface Box :**

Control interface box with process diagram in the front panel and with the same distribution that the different elements located in the unit, for an easy understanding by the student.

All sensors, with their respective signals, are properly manipulated for -10V to +10V computer output.

Sensors connectors in the interface have different pines numbers (from 2 to 16), to avoid connection errors.

Single cable between the control interface box and computer.

The unit control elements are permanently computer controlled, without necessity of changes or connections during the whole process test procedure.

Simultaneously visualization in the computer of all parameters involved in the process.

Calibration of all sensors involved in the process.

Real time curves representation about system responses.

Storage of all the process data and results in a file.

Graphic representation, in real time, of all the process/system responses.

All the actuators' values can be changed at any time from the keyboard allowing the analysis about curves and responses of the whole process.

All the actuators and sensors values and their responses are placed in only one computer screen.

Shield and filtered signals to avoid external interferences.

Real time computer control with flexibility of modifications from the computer keyboard of the parameters, at any moment during the process.

Real time computer control for pumps, compressors, resistances, control valves, etc.

Open control allowing modifications, at any time and in a real time, of parameters involved in the process simultaneously.

Three safety levels, one mechanical in the unit, other electronic in control interface and the third one in the control software.



TRTC/CIB

③ **DAB. Data Acquisition Board:**

PCI Data acquisition board (National Instruments) to be placed in a computer slot. Bus PCI.

Analog input:

Number of channels = 16 single-ended or 8 differential.

Resolution = 16 bits, 1 in 65536.

Sampling rate up to: 250 KS/s (Kilo samples per second).

Input range (V) = ±10V.

Data transfers = DMA, interrupts, programmed I/O. Number of DMA channels = 6.

Analog output:

Number of channels = 2.

Resolution = 16 bits, 1 in 65536.

Maximum output rate up to: 833 KS/s.

Output range (V) = ±10V.

Data transfers = DMA, interrupts, programmed I/O.

Digital Input/Output:

Number of channels = 24 inputs/outputs.

DO or DI Sample Clock frequency: 0 to 1 MHz.

Timing:

Counter/timers = 2.

Resolution: Counter/timers: 32 bits.



DAB

④ **TRTC/CCSOF. Computer Control + Data Acquisition + Data Management Software:**

Compatible with actual Windows operating systems. Graphic and intuitive simulation of the process in screen.

Compatible with the industry standards.

Registration and visualization of all process variables in an automatic and simultaneously way.

Flexible, open and multicontrol software, developed with actual windows graphic systems, acting simultaneously on all process parameters.

Management, processing, comparison and storage of data.

Sampling velocity up to 250,000 data per second guaranteed.

Student calibration system for all sensors involved in the process.

It allows the registration of the alarms state and the graphic representation in real time.

Comparative analysis of the obtained data, after to the process and modification of the conditions during the process.

Open software, allowing to the teacher to modify texts, instructions. Teacher's and student's passwords to facilitate the teacher's control on the student, and allowing the access at different work levels.

This unit allows that the 30 students of the classroom can visualize simultaneously all results and manipulation of the unit, during the process, by using a projector.

⑤ **Cables and Accessories.**

⑥ **Manuals:**

This unit is supplied with 8 manuals: Required Services, Assembly and Installation, Interface and Control Software, Starting-up, Safety, Maintenance, Calibration & Practices Manuals.

* References 1 to 6: TRTC + TRTC/CIB + DAB + TRTC/CCSOF + Cables and Accessories + Manuals are included in the minimum supply, enabling a normal operation.



TRTC/CCSOF



SPECIFICATIONS

Complementary items to the standard supply

PLC. Industrial Control using PLC (7 and 8):

① PLC-PI. PLC Module:

Circuit diagram in the front panel.

Front panel:

Digital inputs (X) and Digital outputs (Y) block:

16 Digital inputs, activated by switches and 16 LEDs for confirmation (red).

14 Digital outputs (through SCSI connector) with 14 LEDs for message (green).

Analog inputs block:

16 Analog inputs (-10V. to + 10V.) (through SCSI connector).

Analog outputs block:

4 Analog outputs (-10V. to + 10V) (through SCSI connector).

Touch screen:

High visibility and multiple functions.

Display of a highly visible status.

Recipe function.

Bar graph function.

Flow display function.

Alarm list.

Multi language function.

True type fonts.

Back panel:

Power supply connector.

Fuse 2A.

RS-232 connector to PC.

Inside:

Power supply outputs: 24 Vdc, 12 Vdc, -12 Vdc, 12 Vdc variable.

Panasonic PLC:

High-speed scan of 0.32 μ sec. for a basic instruction.

Program capacity of 32 Ksteps, with a sufficient comment area.

Free input AC voltage (100 to 240 VAC).

DC input: 16 (24 V DC).

Relay output: 14 (250 VA AC/2 A).

High-speed counter.

Multi-point PID control.

Digital inputs/outputs and analog inputs/outputs Panasonic modules.

Communication RS232 wire, to computer (PC).

② TRTC/PLC-SOF. PLC Control Software:

For this particular unit, always included with PLC supply.



PLC-PI

Items available on request

③ TRTC/CAL. Computer Aided Learning Software (Results Calculation and Analysis).

④ TRTC/FSS. Faults Simulation System.



SACED (EDIBON Computer Control System)

Software Main Screens

Main screen

EDIBON
SENSORS

TEMPERATURE

ST-1	56.47
ST-2	23.16
ST-3	21.83
ST-4	21.32
ST-5	26.64
ST-6	3.94

SW-1 200.22 W
R/min
SR-1 22.61 W/m²
SL-1 0 LUX

ACTUATORS

AR-1 40 60 -80 100

Sensors to Plot

ST-1 ST-2 ST-3 ST-4 ST-5 ST-6 SW-1 SR-1 SL-1 1 GRAPH

Reset Plot Enlarge Plot

AMPLITUDE
TIME (secs)

ST- Temperature sensor, SW- Power sensor, SR- Radiometer, SL- Lux meter, AR- Heating resistance.

Examples of Sensors Calibration screens

ACTUATORS

Analog Input Channel: ST-3
Sensor Name: ST-3
Gain: 107.475 Offset: 4.09457
PTA: 10
Units: 0.002E Calibrated: 4.99E

Simultaneous Calibration

Reference Subject	Sensors	Value	Calibrated	Δ1
<input type="checkbox"/>	ST-1	2.2712	21.5349	23.55
<input checked="" type="checkbox"/>	ST-2	0.3328	25.7898	25.79
<input checked="" type="checkbox"/>	ST-3	0.3328	25.8441	25.84
<input checked="" type="checkbox"/>	ST-4	0.3374	25.5403	25.55
<input checked="" type="checkbox"/>	ST-5	0.3328	25.4378	25.43
<input checked="" type="checkbox"/>	ST-6	0.3408	34.762	34.76
<input type="checkbox"/>	SW-1	0.0000	0.0000	0.00
<input type="checkbox"/>	SR-1	0.0000	0.0000	0.00
<input type="checkbox"/>	SL-1	0.0000	0.0000	0.00
<input type="checkbox"/>	AR-1	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-1	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-2	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-3	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-4	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-5	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-6	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-7	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-8	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-9	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-10	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-11	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-12	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-13	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-14	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-15	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-16	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-17	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-18	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-19	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-20	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-21	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-22	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-23	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-24	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-25	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-26	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-27	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-28	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-29	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-30	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-31	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-32	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-33	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-34	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-35	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-36	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-37	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-38	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-39	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-40	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-41	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-42	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-43	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-44	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-45	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-46	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-47	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-48	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-49	0.0000	0.0000	0.00
<input type="checkbox"/>	AA-50	0.0000	0.0000	0.00

Data table

GAIN	OFFSET	F
ST-1	10.7475	0.0000
ST-2	10.7475	0.0000
ST-3	10.7475	0.0000
ST-4	10.7475	0.0000
ST-5	10.7475	0.0000
ST-6	10.7475	0.0000
SW-1	0.0000	0.0000
SR-1	0.0000	0.0000
SL-1	0.0000	0.0000
AR-1	0.0000	0.0000
AA-1	0.0000	0.0000
AA-2	0.0000	0.0000
AA-3	0.0000	0.0000
AA-4	0.0000	0.0000
AA-5	0.0000	0.0000
AA-6	0.0000	0.0000
AA-7	0.0000	0.0000
AA-8	0.0000	0.0000
AA-9	0.0000	0.0000
AA-10	0.0000	0.0000
AA-11	0.0000	0.0000
AA-12	0.0000	0.0000
AA-13	0.0000	0.0000
AA-14	0.0000	0.0000
AA-15	0.0000	0.0000
AA-16	0.0000	0.0000
AA-17	0.0000	0.0000
AA-18	0.0000	0.0000
AA-19	0.0000	0.0000
AA-20	0.0000	0.0000
AA-21	0.0000	0.0000
AA-22	0.0000	0.0000
AA-23	0.0000	0.0000
AA-24	0.0000	0.0000
AA-25	0.0000	0.0000
AA-26	0.0000	0.0000
AA-27	0.0000	0.0000
AA-28	0.0000	0.0000
AA-29	0.0000	0.0000
AA-30	0.0000	0.0000
AA-31	0.0000	0.0000
AA-32	0.0000	0.0000
AA-33	0.0000	0.0000
AA-34	0.0000	0.0000
AA-35	0.0000	0.0000
AA-36	0.0000	0.0000
AA-37	0.0000	0.0000
AA-38	0.0000	0.0000
AA-39	0.0000	0.0000
AA-40	0.0000	0.0000
AA-41	0.0000	0.0000
AA-42	0.0000	0.0000
AA-43	0.0000	0.0000
AA-44	0.0000	0.0000
AA-45	0.0000	0.0000
AA-46	0.0000	0.0000
AA-47	0.0000	0.0000
AA-48	0.0000	0.0000
AA-49	0.0000	0.0000
AA-50	0.0000	0.0000

ENTER DONE

EXIT SAVE & EXIT




EXERCISES AND PRACTICAL POSSIBILITIES

Practical Possibilities of the Unit:

- 1.- Inverse of the distant square law for the radiation.
- 2.- Stefan-Boltzman 's law.
- 3.- Emission power I.
- 4.- Emission power II.
- 5.- Kirchhoff 's law.
- 6.- Area factors.
- 7.- Inverse of the distant square law for the light.
- 8.- A ' lambert cosine 's law.
- 9.- A ' lambert absorption 's law.
- Other possible practices:
- 10.-Sensors Calibration.
- Practices to be done by PLC Module (PLC-Pi) + PLC Control Software:
- 11.-Control of the TRTC unit process through the control interface box without the computer.
- 12.-Visualization of all the sensors values used in the TRTC unit process.
- 13.-Calibration of all sensors included in the TRTC unit process.
- 14.-Hand on of all the actuators involved in the TRTC unit process.
- 15.-Realization of different experiments, in automatic way, without having in front the unit. [This experiment can be decided previously].
- 16.- Simulation of outside actions, in the cases do not exist hardware elements.(Example: test of complementary tanks, complementary industrial environment to the process to be studied, etc).
- 17.- PLC hardware general use and manipulation.
- 18.- PLC process application for TRTC unit.
- 19.- PLC astructure.
- 20.- PLC inputs and outputs configuration.
- 21.- PLC configuration possibilities.
- 22.- PLC program languages.
- 23.- PLC different programming standard languages .
- 24.- New configuration and development of new process.
- 25.- Hand on an established process.
- 26.- To visualize and see the results and to make comparisons with the TRTC unit process.
- 27.- Possibility of creating new process in relation with the TRTC unit.
- 28.- PLC Programming Exercises.
- 29.- Own PLC applications in accordance with teacher and student requirements.

POSSIBILITIES OF OTHER AVAILABLE EXPANSIONS



ORDER INFORMATION

Items supplied as standard

Minimum configuration for normal operation includes:

- ① Unit: TRTC, Thermal Radiation Unit.
- ② TRTC/CIB, Control Interface Box.
- ③ DAB, Data Acquisition Board.
- ④ TRTC/CCSOF, Computer Control + Data Acquisition + Data Management Software.
- ⑤ Cables and Accessories.
- ⑥ Manuals.

* **IMPORTANT:** Under TRTC we always supply all the elements for immediate running as 1, 2, 3, 4, 5 and 6.

Complementary items to the standard supply

- PLC, Industrial Control using PLC (7 and 8):
- ⑦ PCL-PLC Module.
 - ⑧ TRTC/PLC-SOF, PLC Control Software.
 - ⑨ TRTC/CAL, Computer Aided Learning Software (Results Calculation and Analysis). (Available on request).
 - ⑩ TRTC/FSS, Faults Simulation System. (Available on request).
- Expansions**
- ⑪ Mini ESN, Multipost EDIBON Mini Scada-Net System.
 - ⑫ ESN, Multipost EDIBON Scada-Net System.

REQUIRED SERVICES

- Electrical supply: 220V/50 Hz or 110V/60 Hz.
- Computer (PC).

DIMENSIONS & WEIGHTS

- TRTC Unit: -Dimensions: 1400 x 820 x 900 mm. approx.
-Weight: 40 Kg. approx.
- Control Interface Box: -Dimensions: 490 x 330 x 310 mm. approx.
-Weight: 10 Kg. approx.
- PLC Module (PLC-Pi): -Dimensions: 490 x 330 x 310 mm. approx.
-Weight: 30 Kg. approx.

* Specifications subject to change without previous notice, due to the convenience of improvements of the product.

EDIBON International, S.A.

REPRESENTATIVE: